

EXPERIMENT NO. 7

Aim: To determine the Total Hardness in given water sample by complexometric titration method.

Apparatus required: Burette, conical flask, pipette, beaker and Funnel.

Chemical required: 0.1 M- EDTA, 0.01 M - CaCO_3 , EBT Indicator ,buffer solution (NH_3 / NH_4Cl) and given water sample.

Theory: The hardness of water can be determined by complexometric titration and EDTA is used as complexing agents. The Ca^{2+} and Mg^{2+} are present in water and titrated with EDTA using as a EBT indicator.

Chemical Equations:

- I) The EBT indicator forms an unstable complex with Ca^{2+} or Mg^{2+} at pH of 9 to 10. The complex is wine red in colour.
- $$\text{Ca}^{2+} / \text{Mg}^{2+} + \text{EBT} \longrightarrow [\text{Ca} / \text{Mg} - \text{EBT}]$$
- (Unstable complex – Wine red)
- II) As this solution is titrated against EDTA, the free Ca^{2+} or Mg^{2+} ions in water forms stable metal ion EDTA complex.
- $$\text{Ca}^{2+} / \text{Mg}^{2+} + \text{EDTA} \longrightarrow [\text{Ca} / \text{Mg} - \text{EDTA}]$$
- (Stable complex - Colourless)
- III) The free metal ions are complexed, the EDTA replaces Ca or Mg ions from the unstable indicator and the indicator is free. Since the indicator is blue. So, the end point is appearance of blue colour.
- $$[\text{Ca} / \text{Mg} - \text{EBT}] + \text{EDTA} \rightarrow [\text{Ca} / \text{Mg} - \text{EDTA}] + \text{EBT}$$
- Wine red Colourless Blue

Procedure :

i) Standardization of EDTA / CaCO_3

1. EDTA is filled in burette.
2. 10 ml CaCO_3 is pipette out in a conical flask.
3. 5 ml of buffer is added in the conical flask.
4. 2-3 drops of EBT (Indicator) is also added.
5. Titrate against the EDTA solution till the colour of the solution changes from wine red to blue.
6. Note the final reading and repeat to get three concordant readings.

ii) EDTA Vs Water Sample

1. EDTA solution is filled in the burette.
2. 10 ml of water sample is pipette out in conical flask.
3. Add 5 ml of buffer solution and 2-3 drops of EBT indicator.
4. Titrate against the EDTA solution till the colour of the solution changes from wine red to blue.
5. Write the final reading and repeat to get three concordant readings.

Observations:

i) Standardization of EDTA

Solution in burette = EDTA
In conical flask = CaCO_3 (10 ml)
Indicator = EBT
End Point = Wine Red to Blue

Observation Table:

Serial no.	Burette Reading		Volume of EDTA (ml)
	Initial	Final	
1	0	11.5	11.5
2	11.5	23	11.5

ii) EDTA Vs Water Sample

Solution in burette = EDTA
Solution in flask = Water Sample (10 ml)
Indicator = EBT
End Point = Wine Red to Blue

Observation Table:

Serial no.	Burette Reading		Volume of EDTA (ml)
	Initial	Final	
1	0	12.5	12.5
2	12.5	25	12.5

Calculations :

i) Standardization of EDTA (EDTA/CaCO₃)

$$\begin{aligned}M_1 V_1 &= M_2 V_2 \\(\text{EDTA}) &(\text{CaCO}_3) \\M_1 &= M_2 V_2 / V_1 \\M_{1(\text{EDTA})} &= 0.1 \times 10 / 11.5 \\&= 0.086 \text{ M}\end{aligned}$$

iii) EDTA Vs Water sample

$$\begin{aligned}M_1 V_1 &= M_2 V_2 \\(\text{Water}) &(\text{EDTA}) \\M_{1(\text{water})} &= M_2 V_2 / V_1 \\&= 0.086 \times 12.5 / 10 \\&= 0.107 \text{ M} \\ \text{Strength} &= \text{Molarity of hard water} \times \text{Molecular mass of CaCO}_3 \\&= 0.107 \times 100 \text{ g / L} \\&= 10.7 \times 1000 \\&= 1070 \text{ ppm}\end{aligned}$$


Result: Total Hardness of water sample is 1070 ppm.

Precautions:

1. The burette, pipette and conical flask should be rinsed with distilled water.
2. The colour change near the end point is very slow and thus should be observed carefully.

Structure:

EBT (Eriochrome Black T) Molecular formula- C₂₀H₁₂N₃NaO₇S



EDTA (Ethylene Diamine Tetra Acitic Acid)

Chemical formula: $\text{C}_{10}\text{H}_{16}\text{N}_2\text{O}_8$

